

Ring Test Results

Test	Assigned value (µg/kg)	Result (µg/kg)	Z-score
A.I.A. Ring Test Lotto RTB1 240919 Aflatoxin B1 in Corn September 2019	<limit of detection	<limit of detection	-0.142
	7.2	7.15	
	13.9	12.90	
	22.3	22.07	
FAPAS Food Chemistry Proficiency Test 04375 Mycotoxin Contamination in Animal Feed (Cereal Based) December 2019 - January 2020	14.6	12.78	-0.6
FAPAS Food Chemistry Proficiency Test 04378 Multi-Mycotoxins in Maize Flour February - March 2020	6.65	6.48	-0.1
FAPAS Food Chemistry Proficiency Test 04384 Multi-Mycotoxins in Maize Flour April - May 2020	6.13	5.95	-0.1
FAPAS Food Chemistry Proficiency Test 04390 Mycotoxins in Hazelnut (Water/Nut Slurry) July - August 2020	4.50	3.67	-0.8



FAPAS Food Chemistry Proficiency Test 04402 Mycotoxin Contamination in Animal Feed (Cereal Based) December 2020 – January 2021	20.4	18.18	-0.5
FAPAS Food Chemistry Proficiency Test 04408 Aflatoxins in Peanut (Water/Nut Slurry) March – April 2021	4.86	5.64	0.7
FAPAS Food Chemistry Proficiency Test 04411 Multi-Mycotoxins in Maize Flour April – June 2021	5.35	5.25	-0.1
FAPAS Food Chemistry Proficiency Test 04417 Mycotoxins in Hazelnut (Water/Nut Slurry) June - August 2021	4.20	4.81	0.7
FAPAS Food Chemistry Proficiency Test 04418 Multi-Mycotoxins in Maize Flour July – September 2021	3.12	3.35	0.3

