



BIO-SHIELD FUMONISIN 5

Cat.Number B5548/B5596

Ring Test Results

Test	Assigned value ($\mu\text{g}/\text{kg}$)	Result ($\mu\text{g}/\text{kg}$)	Z-score
FAPAS Food Chemistry Proficiency Test 04384 Multi-Mycotoxins in Maize Flour April - May 2020	752	612.19	-1.1
FAPAS Food Chemistry Proficiency Test 04411 Multi-Mycotoxins in Maize Flour April - June 2021	962	1015.13	0.3
FAPAS Food Chemistry Proficiency Test 04418 Multi - Mycotoxins in Maize Flour July - September 2021	720	658.5	-0.5